

# DIRECT DRUG TESTING LTD



Contact: Jo Kirk • Telephone: 027 243 2544 • E-mail: jo@directdrugtesting.co.nz

## ProScreen™ 6 Panel Dip with Adulterants

New Zealand's **Best Selling**  
Drug Testing Dip Test!



| ADULTERATION COLOR COMPARISON CHART   |                |                     |  |                 |
|---------------------------------------|----------------|---------------------|--|-----------------|
| TEST AND READING TIME                 | ABNORMAL (LOW) | NORMAL              |  | ABNORMAL (HIGH) |
| Creatinine (Cr)<br>45 seconds         | <10 10         | 20 50 100 200 mg/dl |  |                 |
| Nitrite (Ni)<br>45 seconds            |                | 0 0.1-0.2 0.5-0.8   |  | > 15 mg/dl      |
| pH<br>Immediate                       | 2 3            | 4 5 7 9             |  |                 |
| Bleach (BI)<br>45 seconds             |                | Negative            |  | Positive        |
| Specific Gravity (S.G.)<br>60 seconds | 1.000          | 1.005 1.015 1.025   |  | ≥ 1.030         |



### TEST FOR 6 DRUGS

Built in Specimen Validity test to defeat specimen adulteration.

Eliminates Collector Error – the only no flood drug test!

Fast.  
Results typically ready in 1-2 minutes.

Accurate.  
Verified to the Australian Standards and used by laboratories, government agencies and corporations across New Zealand!

Affordable.  
Best priced and the most features of any dip test in New Zealand — Guaranteed!

### Drug Cutoff Levels

|      |           |
|------|-----------|
| THC  | 50 ng/ml  |
| COC  | 300 ng/ml |
| AMP  | 300 ng/ml |
| METH | 300 ng/ml |
| OPI  | 300 ng/ml |
| BZO  | 200 ng/ml |